

*John Sterling, Inc.*


**“Strength through Knowledge”**

## **New Product for the Medical Gas Market**

### **Medi-Purge™ a patent pending NF nitrogen purge device**


**(part #MP-A)**

- Easily adapts to Chemetron® medical gas wall or console outlets / inlets
- Screws into outlet / inlet back body
- Allows NF nitrogen into an outlet back body
- Allows NF nitrogen to flow out of an outlet back body
- Material - brass - oxygen cleaned and bagged
- 1/4” hole for O2 analyzer to ensure all air has been expelled
- 1/4” male thread for NF nitrogen line
- Flow allowances of 37.2 SCFM @ 35 psig
- Easy way to fill lines with NF nitrogen for pipe that won't be worked on for more than 8 hours
- Braze multiple outlets quickly by placing a Medi-Purge™ device on all same gas outlets in a zone



Screwing the threaded plunger allows easy access to flow NF nitrogen in or out of the outlet. (The Medi-Purge™ will open check on console outlets)

Through hole to allow for NF nitrogen to flow into or out of an outlet / inlet



1/4” hole for oxygen analyzer

1/4” male NPT threads for NF Nitrogen connection

Fits console or wall outlets

[www.medi-purge.com](http://www.medi-purge.com)

Ph: (859) 586-2971

Fax: (859) 586-4207

2940 Hebron Park Dr. Suite 201 Hebron, KY 41048